



qLabs[®] Q3 PRO ELECTROMETER

for PT-INR/APTT/FIB results on whole blood

“Treatment decisions
are often based on
**Samples that are
hours old**”
-ICU Physician



Stay Confident Using **Q3 PRO Electrometer** for near patient testing

OVERVIEW:

Q3 PRO electrometer is suitable for patients as a self-testing device who are on OAT treatment which is essential in overall prophylaxis and treatment for thrombotic disorders.

Q3 PRO also finds its use across:

- | | | | |
|----------------------|------------|-------|--------------------|
| • Stroke/INR clinics | Small labs | ICU's | GPO's |
| • Cardiac OPD's | Casualty | ITU's | Hemophilia clinics |

TECHNICAL SPECIFICATION

- **Display** : 3.5" Color LCD
- **Sample Type** : Fresh whole blood (Capillary & Venous)
- **Sample Volume** : Less than 10µl
- **Temperature** : Inbuilt temperature
- **Quality** : Build-in quality control
- **Language** : Multi-language menu
- **Measurable Range** : INR 0.5 to 8.0
- **Memory Capacity** : Up to 2000 records
- **Dimensions** : 148 x 70 x 26 (mm)
- **Weight** : 199 g
- **Temperature** : 10 °C - 35 °C
- **Relative humidity** : 10%- 90%
- **Power Supply** : 100V-240V AC
- **Battery** : 2500mAH lithium built-in battery
- **Backup** : Up-to 120 hours of usages



CURRENT TEST MENU

- PT-INR/APTT/FIB



Follow Us On

- www.operonbiotech.com
 info@operonbiotech.com
 @operonbio1 @operonbio1
 @operonbio1 @operonbio1

Operon Biotech & Healthcare

- GVR Hexa, No. 428, 3rd Floor, 8th Main,
Koramangala - 4th Block, Bangalore -
560047
 080-42027704